

Manufacturer	SAFEMED Gil, Bluszcz Spółka Jawna
Contact	Ul. Wielicka 181A/22, 30-663 Kraków Tel. +48 12 681 42 52 fax +48 12 398 41 52 email: biuro@safemed.pl



EFFECT PF BLUE

nitrile examination and protective gloves, powder-free, disposable, non-sterile	
Product code	R-N-B-EF100-size
Material	Nitrile butadiene rubber
External surface	Textured fingers
Internal surface	Powder-free, chlorinated
Colour	Blue
Shape	Ambidextrous
Cuff	Rolled
AQL	1.5 according to G1 i ISO 2859
Classification	Medical device, class I in accordance with MDR 2017/745, Personal Protective Equipment category III in accordance with Regulation (EU) 2016/425
Standards	EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 1041, EN ISO 14971, EN ISO 15223-1, EN 420, EN ISO 374-1, EN ISO 374-5, EN 374-2, EN 374-4, EN 16523-1
Food contact	Approved, in accordance with Regulation 10/2011, 1935/2004
Quality control system	ISO 9001, ISO 13485
Packaging	a'100/ 10 x a'100
Storage conditions	Temp. 5°C – 38°C
Sizes	XS, S, M, L, XL

PHYSICAL PROPERTIES

Single wall thickness	Finger: min. 0.07 mm, palm: min. 0.06 mm, cuff: min. 0.04 mm
Length	Min. 240 mm
Powder content	< 2 mg/glove
Force at break before/ and after aging:	min. 6N / min. 6N

Date of issue: 27.05.2021